

1/24/2017

Installing Xcalibur 4.0 with Foundation 3.1 SP2 on Windows 10 pro

For reviewing data on a Personnel Workstation you need to manually install the MSMQ feature from control panel, and then install the foundation and Xcalibur. (Note, this is not intended for instrument PC's, Win 10 is not yet supported by Thermo.)

http://portal.thermo-brims.com/index.php?option=com_discussions&view=thread&catid=4%3Aproteomics-forum&thread=3242%3Axcalibur-4-0-error-mqrt-dll&Itemid=54

Installing Message Queuing (MSMQ)

1. Open **Control Panel**.
2. Click **Programs** and then, under **Programs and Features**, click **Turn Windows Features on and off**.
3. Expand Microsoft Message Queue (MSMQ) Server, expand Microsoft Message Queue (MSMQ) Server Core, and then select the check boxes for the following Message Queuing features to install:
 - MSMQ Active Directory Domain Services Integration (for computers joined to a Domain).
 - MSMQ HTTP Support
4. Click **OK**.
5. If you are prompted to restart the computer, click **OK** to complete the installation. (I restarted it anyway, w/o prompt)

Then install Xcalibur 4.0: it will automatically install foundation, reboot, and then install Xcalibur and Freestyle

Seems to work:

The image contains three screenshots. The leftmost screenshot shows the Thermo Xcalibur Qual Browser interface with two mass spectra plots. The top plot shows relative abundance versus time (min) with peaks labeled at 5.87, 7.57, 26.50, 38.89, 41.86, 49.73, 53.10, 55.49, 61.33, 70.20, 73.58, 74.56, 78.46, 90.38, 95.40, 96.96, 105.72, and 113.21. The bottom plot shows relative abundance versus m/z with peaks labeled at 415.07, 476.91, 555.64, 658.41, 727.25, 813.46, 889.39, 905.46, 993.06, and 1093.30. The middle screenshot shows the Windows 10 Control Panel 'System' window, displaying basic information about the computer: Windows edition (Windows 10 Pro), processor (Intel(R) Core(TM) i7-6700K CPU @ 4.00GHz), installed memory (24.0 GB), and system type (64-bit Operating System). The rightmost screenshot shows the 'Version Info' window under 'Installed Software', listing the following installed software: Foundation 3.1 SP2 (Release Date: November 18, 2015; Install Date: January 24, 2017; Description: Thermo Foundation 3.1 Release; Version: 3.1.94.0), Freestyle 1.1 SP1 (Release Date: March 3, 2015; Install Date: January 24, 2017; Description: Thermo Freestyle 1.1 Release; Build Number: 175; Revision Number: 0), Xcalibur 4.0 (Release Date: November 18, 2015; Install Date: January 24, 2017; Description: Thermo Xcalibur Release; Version: 4.0.27.13), and Configured Instruments (No instruments configured).